

| MODEL | CODE | CAPACITY | DIVISION | WEIGHT SURFACE |
|-------|-------|----------|----------|----------------|
| C201 | 00008 | 200kg | 100g | 275 x 530mm |

Accessories Included: Telescopic height rod measuring range 60cm~212cm, Adjustable feet, Wheels to move; Warranty 1 year


| SERVICES: | CODE |
|---------------------------|--------|
| Verification | AOMOL |
| Calibration's certificate | ACOSIT |
| Warranty extension | AGAES |

FUNCTIONS

Weight, height rod measuring, Zero adjustable

PRODUCT DETAILS

- Certifications:** EU directive NAWI approved class III.
EEC directive MDD risk class Im, with measurement function.
- Structure:** Metal painted structure with rubber platform. Weight system with lever mechanism with knife-edges. Weight easy and clear to read.
- Dimensions mm:** Total (w)450x(L)600(h)1500,
Platform (w)275x(L)530x(h)115, Pole (h)1200.
- Packaging:** 1pcs dimension (w)525 x(L)1500x(h)170mm,
Net weight 15kg, Gross Weight 17.5kg



Wheels to move



Telescopic height rod



Lock beam scale